



wave•one™

Clinical Procedure Card



- The wave•one™ instruments must only be used in reciprocating motion with appropriate Dentsply Maillefer wave•one™ motors
- For more information please refer to the wave•one™ DFU





wave•one™

Reciprocating files



wave•one™ **small** reciprocating file

REF A 0700 2210 0 S01- **3x** small (021 .06) - 21mm + DFU

REF A 0700 2250 0 S01- **3x** small (021 .06) - 25mm + DFU

REF A 0700 2310 0 S01- **3x** small (021 .06) - 31mm + DFU



wave•one™ **primary** reciprocating file

REF A 0700 2210 0 P01- **3x** primary (025 .08) - 21mm + DFU

REF A 0700 2250 0 P01- **3x** primary (025 .08) - 25mm + DFU

REF A 0700 2310 0 P01- **3x** primary (025 .08) - 31mm + DFU



wave•one™ **large** reciprocating file

REF A 0700 2210 0 L01- **3x** large (040 .08) - 21mm + DFU

REF A 0700 2250 0 L01- **3x** large (040 .08) - 25mm + DFU

REF A 0700 2310 0 L01- **3x** large (040 .08) - 31mm + DFU



wave•one™ **assortment** reciprocating file

REF A 0700 2219 0A01 assorted - 21mm

1x small (021 .06) + **1x** primary (025 .08) + **1x** large (040 .08) + DFU

REF A 0700 2259 0A01 assorted - 25mm

1x small (021 .06) + **1x** primary (025 .08) + **1x** large (040 .08) + DFU

REF A 0700 2319 0A01 assorted - 31mm

1x small (021 .06) + **1x** primary (025 .08) + **1x** large (040 .08) + DFU

- Only one NiTi instrument per root canal in most cases
- Decreases the global shaping time by up to 40%
- Reduced risk of a screwing effect and file breakage*
- Flexibility and greater resistance to cyclic fatigue, the leading cause of file separation
- Respects the root canal anatomy
- Single use as new standard of care



*Compared to leading brands

+
**WE
KNOW
ENDO.**

CE
0459



Dentsply Maillefer • Chemin du Verger 3 • CH-1338 Ballaigues • Switzerland • www.dentsplymaillefer.com